


*Volumup* *Lifting*  
**V.V.L**



## What is PDO Thread

PDO (Polydioxanone) is a synthetic and biodegradable monofilament suture. Due to its safety and effectiveness, it has been used widely. Inserted in to skin, PDO stimulates fibroblast and induce micro circulation.





## Volumup Lifting Effect

- Lifting
- Skin Rejuvenation (Brightness & elasticity)
- Wrinkle Removal



# Volumup *V*Lifting

**Mono**



Gauge	Needle length	Thread length	USP
30G	30mm	40mm	6-0
29G	40mm	60mm	6-0
29G	50mm	80mm	6-0
27G	60mm	90mm	5-0

※ More sizes are available.




# Volumup V Lifting **Screw**



Gauge	Needle length	Thread length	USP
30G	30mm	40mm	6-0
30G	40mm	60mm	6-0
27G	60mm	90mm	5-0

※ More sizes are available.

# Volumup V Lifting **Cog**



Bi-directional cog holds skin tightly for lifting.




Gauge	Needle length	Thread length	USP
21G	60mm	130mm	3-0
21G	90mm	150mm	3-0
23G	60mm	130mm	3-0
23G	90mm	150mm	3-0


※ More sizes are available.



Adoption of blunt cannula  
minimize any possible  
pain & side effects dramatically.



Blunt cannula guarantees  
safe & easy treatment.



Gauge	Needle length	Thread length	USP
21G	60mm	130mm	3-0
21G	90mm	150mm	3-0
23G	60mm	130mm	3-0
23G	90mm	150mm	3-0

※ More size are available.





# Volumup Lifting

## Canula Mono

Product Description Micro cannula mono with blunt tip

Specially designed for eye zone treatment

Blunt tip cannula type needle provides safe, easy and effective eye zone treatment

Reduce downtime, bruising, pain and complaints

Less ecchymosis during treatment for eye zone

Less bruising, ecchymosis, and pain with faster recovery



Gauge	Needle length	Thread length
30G	28,5mm	40mm

# Volumup Lifting

## Super Screw

Major features of Super Screw

Elastic force effect by screw shape PDO

Continuous physical stimulation to skin

Increase composition of collagen

Increase secretion of Growth factors and cytokines

High density linear effect by fibrous bridging effect

Augmentation effect by insertion of densely coiled thread (coiled length 2cm, spread length 10cm)

Continuous mechanotransduction effect

High cell cultivation with large contact area



Gauge	Needle length	Thread length
28G	40mm	100mm






## PLLA

(Mono / Cog type available)



PLLA is a safe and absorbable material which has been used for wrinkle treatment fillers for a long time. PLLA suture stimulates the formation of collagen continuously in the skin to recover your facelift and improve wrinkles. It last about 18–24 months before hydrolysis process.



## PCL

(Only Mono type available)

PCL suture helps our skin produce collagen, which is responsible for skin resilience and regeneration due to its excellent characteristics. PCL suture is widely used as surgery stitch thread which is also safe to human body. It is acknowledged to be safe product approved by FDA. PCL suture sets you free from frequent treatments as it lasts two years. PCL is safe because it is completely decomposed and absorbed in our body through hydrolysis process.



# Why is the *Volumup* Lifting

## Long lasting Quality

According to PDOs' characteristics of water absorption, most other companies products' shelf life time is only one year.

**Volumup V lifting is using dry holds unique vaccum dry technology.**

We guarantees 2 years life time.

It leads to more effective and safe treatment.

## Originality

PDO thread was invented in 1995 Dr. Kim Dong Jun in our company by and has been developing technology since.

Holder of patents in Korea and China.

## Certificate

- CE Certificate for models including cog type
- ISO13485
- Patent No. 10-0803388
- Patent No. 200710161694.7

